# Novel Achiral Building Blocks and Solid Supports for Preparation of Multi-Labeled Oligonucleotides



Khirud Gogoi, Vladimir Yu. Vvedensky and Andrei P. Guzaev

AM Chemicals LLC, 4065 Oceanside Blvd., Ste M, Oceanside, CA, 92056, USA khirudg@gmail.com, aguzaev@amchemicals.com



#### Abstract

Novel achiral non-nucleosidic building blocks for the efficient preparation of single- or multi- labeled oligonucleotides have been developed. Incorporation of a wide variety of ligands including GalNAc cluster, fluorescent labels and quenchers, biotin, hydrophobic moieties, azido and alkyne groups suitable for the use in click-type conjugation has been demonstrated.

Synthesis of Azetidine- and Piperidine core structures and attachement of various ligands



#### Synthesis & deprotection of oligonucleotides

Sequence 1: 5'-TTT TTT TTT TTT TTT TTT TTT TTT Ligand 3' Sequence 2: 5'-TAG TGC TAG ATG CCT-Ligand 3' Sequence 3: 5'-CCA CTA CCT GAG CAC CCA GTT-Ligand 3' Sequence 4: 5'-CTG GGT GCT CAG GTA GTG GTT-Ligand 3' Sequence 5: 5'-Ligand-CCA CTA CCT GAG CAC CCA GTT 3' Sequence 6: 5'-Ligand-CCA CTA CCT GAG CAC CCA GTT-Ligand 3' Sequence 7: 5'-CsCsAs CsTsAs CsCsTs GsAsGs CsAsCs CsCsAs GsTsT-Ligand 3' Sequence 8: 5'-r(CCA CUA CCU GAG CAC CCA GUU)-Ligand 3' Sequence 9: 5'-r(CsCsAs CsUsAs CsCsUs GsAsGs CsAsCs CsCsAs GsUsU)-Ligand 3'

## Introduction

- Non-nucleosidic building blocks used in oligonucleotide chemistry for non-radioactive labeling and for attachment of various ligands should result in the attachment of a linker stable under the conditions of oligonucleotide chain assembly and deprotection by most common basic deprotection agents.
- The currently available building blocks derived from 1,2-diols suffer an unwanted loss of the linker due to cleavage of phosphodiester group.
- iii) Generally the side chain of a building block is attached to the core structure via an amido function. When the amido group is placed in proximity to the nearest phosphodiester moiety, a minor loss of the entire linker via an intramolecular cyclization may occur and an oligonucleotide bearing a terminal phosphate group is formed
- iv) The attached linker should not generate new chiral centers in the oligonucleotide conjugate obtained. Ideally, the building block itself should not possess any chiral or prochiral carbon atoms.



Base-catalyzed dephosphorylation is always ii) particularly rapid when the conformation is iii) particularly rapid when phosphate is



#### Synthesis of Solid Supports and Phosphoramidite Building Blocks



**15**: n = 2;

**16**: n = 1

) Piperidine, MeOH

ii) D-Biotin, TBTU, iPr<sub>2</sub>NET

#### <u>Cleavage and deprotection of the oligonucleotides:</u>

i) Aqueous NH<sub>3</sub> (28%), 60<sup>o</sup>C, 8-12 hour.

ii) Acetonitrile-Diethylamine (5:1) treatment for 3 minutes, followed by treatment with aqueous NH<sub>3</sub> (28%), 60°C, 8-12 hour.

iii) 1:1 AMA (Aqueous NH<sub>3</sub>: Methylamine), 60°C, 15 minutes

iv) Acetonitrile-Diethylamine (5:1) treatment for 3 minutes, followed by treatment with Toluene-Ethylenediamine (1:1) for 2 hours at room temperature.

v) 50 mM K<sub>2</sub>CO<sub>3</sub> in Methanol at room temperature (45 °C for CF<sub>3</sub>CO protected Amine oligonucleotides)

vi) Acetonitrile-Diethylamine (5:1) treatment for 3 minutes, followed by treatment with 50 mM K<sub>2</sub>CO<sub>3</sub> in Methanol at room temperature (45 °C for CF<sub>3</sub>CO protected Amine oligonucleotides)

## Results:

All the Methods (i) to (vi) were suitable for oligonucleotides of Alkyne, Dabsyl, Dabcyl, Pyrene, Fluorescein, Cholesterol and GalNAc.

Method (ii) to (vi ) were suitable for deprotection of Biotin labeled oligonucleotides. Method (i), (ii) & (vi, 45 °C) lead to complete deprotection of the CF<sub>3</sub>CO protecting group of the Amine labeled oligonucleotides.

#### HPLC profiles of Ligand-oligonucleotides conjugates



